

PROCESSING OF MACROALGAL SAMPLES

TRAITEMENT DES ÉCHANTILLONS DE MACRO-ALGUES



LIST OF EQUIPMENTS

List of equipments
Liste du matériel nécessaire



- 125 µm, 40 µm and 20 µm mesh size sieves
- 50mL screw cap tubes and rack
- A thin indelible marker
- A reservoir filled with seawater collected *in situ* from sampling sites
- A squeeze bottle with filtered seawater
- A small funnel
- x2 100mL graduated cylinders
- Paraffin tape
- 5mL plastic pipettes and pump

List of equipments
Liste du matériel nécessaire



- A notebook and a pen
- Formalin solution - 37% weight in water
- A micropipet (P1000) with appropriate disposable tips
- Plastic Pasteur pipettes
- A portable scale
- A brush 3cm wide