

PROCESSING OF WINDOW SCREEN SAMPLES

TRAITEMENT DES ÉCHANTILLONS DE WINDOW SCREEN



LIST OF EQUIPMENTS

List of equipments
Liste du matériel nécessaire



- A 1L graduated cylinder
- A 200-300 μm pore size mesh sieve
- A funnel large enough to fit the 200-300 μm sieve
- Some 20 μm nylon mesh pieces cut into 4 cm squares
- A filtration unit composed of a 200mL or 500mL filter funnel, a 25 mm diameter filter base and its stopper, and a 500mL or 1L side-arm flask
- A hand vacuum pump
- 2 pairs of flat forceps
- A squeeze bottle with filtered seawater

List of equipments
Liste du matériel nécessaire



- 15mL conical screw cap tubes and rack
- Formalin solution - 37% weight in water
- Pure ethanol
- 10mL plastic pipettes in individual packaging
- A micropipet (P1000) with appropriate disposable tips
- Plastic Pasteur pipettes
- A thin indelible marker
- A notebook and a pen
- Paraffin tape